

How to Use Nasal Atomizer



Step 1: Arrange materials from package. **Carefully read instructions** on prescribed medication label before continuing.


Disclaimer: Your package may not contain some products appearing in this picture.

Step 2: Add 10mL of medication to medicine cup and fill the syringe.



Step 3: Attach the nasal atomizer tip to the syringe by screwing on.

Step 4. Place the tip of the atomizer in the nostril. Press the plunger and gently spray half (5mL) of the medication into the nostril. Move the device over to the opposite nostril and repeat.



Cleaning:

1. Clean the device with dishwashing detergent and hot water. After cleaning, rinse thoroughly with distilled water and allow to air dry.
2. Allow to dry between each use.
3. If clogged, submerge tip in warm water and pull back on plunger to release clogged particles. Remove tip and empty syringe.

Refills: Don't forget, you may have refills! Simply text your name and prescription number to (910) 763-1896 to refill your prescription. Or ask us about auto-refills!

*Please fill out our Patient Survey found on our website under the contact menu.
We value your feedback.*

Medical Center Specialty Pharmacy
912 S. 16th Street
Wilmington, NC 28401
(910)763-1896 / (800) 930-8785



www.mcspnc.com