

PATIENT INFORMATION SHEET

HOW TO USE A NASONEB

INSTRUCTIONS:



1. Arrange materials from the package. **Carefully read the instructions** on prescribed medication labeling before continuing. Disclaimer: Your package may not contain some products appearing in this picture.
2. Place half of the specified dose of medication in the NasoNeb® cup as directed on prescription bottle.
3. Tilt your head slightly downward and place the NasoNeb nozzle into one nostril, aim it straight in or slightly away from the septum. Turn the compressor on, and turn off once the cup is empty
4. Repeat **Step 3** in your other nostril until the cup is empty or 1 minute has passed. Aim the NasoNeb Cup straight in your nostril or slightly away from the septum. **Do not** aim it directly at the septum which divides the two nostrils.

**Note: Do no use Nasoneb Moisturizing Solution with medication.*



TECHNIQUE: TILT YOUR HEAD SLIGHTLY DOWNWARD. BREATHE THROUGH BOTH YOUR NOSE AND MOUTH WITH THE COMPRESSOR ON. THE TIP OF THE INSERT SHOULD BARELY PENETRATE YOUR NOSTRIL. TURN OFF COMPRESSOR BEFORE REMOVING THE MEDICINE CUP FROM YOUR NOSTRIL. SOME MEDICATIONS MAY LEAVE FOAM IN THE CUP. THIS IS EXPECTED.

CLEANING

Clean the nasal nebulizer with warm, soapy water after each use. Avoid the use of harsh detergents or cleaning agents. Ensure the compressor is turned off and unplugged and the product is disassembled.

For concerns related to the function of the Nasonebulizer, reach out to the manufacturer at: 800-833-9653. Their team is diligent in helping troubleshoot issues with the machine as they arise.

Address:

Medical Center Pharmacy
912 S. 16th Street Wilmington, NC
28401

Contact:

(910)763-1896 / (800) 930-8785

Follow us on social media:



@mcspwilmington

